

HEALTH IS THE HEART OF OUR BUSINESS.

HOLTER-RECORDER EP820 / EP820-12


DIGITAL HOLTER RECORDER WITH PM DETECTION

High Patient Acceptance

3 or 12 Recording Channels

Recording Time up to 7 Days

Real Time ECG Display



MANAGE COMPLEX DATA EASILY

AMEDTEC ECG^{pro}

TECHNICAL DATA

Dimensions	Size	68 mm x 53 mm x 16 mm (L x W x H)
	Weight	42 g
	Operating Position	any orientation
Electrical	Gain Settings	0,5 x, 1 x and 2 x
	Connector	26 PIN
	Patient Cable	4, 5 or 10 electrodes
Functional	Recording Channels	3 or 12
	Sensitivity	2,5 µV/ LSB
	Recording	full disclosure
	Data Transfer	USB 2.0 or Card Reader
	Sample Rate	10.000/s (PM-Detection)
	Frequency Response	0,05 Hz to 60 Hz, -3 dB
	ECG Signal View	LCD at Hook-Up or on Demand
	PM Detection	permanently
Battery	Battery Type	1 AAA Battery, 1,5 V Alkaline, Lithium or NiMH
	Life (3 Channels)	Alkaline: ≥ 96 h Lithium: ≥ 168 h
	Life (12 Channels)	Alkaline: ≥ 48 h


ACCESSORIES



Holter ECG electrodes



disposable pouch




disposable sanitary
neck strap



patient cable

AMEDTEC Medizintechnik Aue GmbH
Schneeberger Straße 5 · 08280 Aue · Germany
Phone: +49 (0)3771 59 82 70
Fax: +49 (0)3771 59 82 790
info@amedtec.de · www.amedtec.de



CE Certificate:
0494 DIN EN ISO 13485

Errors and omissions excepted
M0672 019-0842-2_02