

HEALTH IS THE HEART OF OUR BUSINESS.

QUICKELS DECAPUS® III ECG VACUUM SYSTEM



A SYSTEM CARING FOR BEST POSSIBLE HYGIENIC SAFETY AND HIGHEST EFFICIENCY AT ALL TIMES USING **THE CLICK-IN AND CLICK-OUT FEATURE**

CLICK-IN / CLICK-OUT - Efficient and economical to use by saving time when applying

Highest hygiene safety standard while using disposable electrodes

High lead quality providing reliable and accurate ECG diagnostics

Mobil and easy to use



MANAGE COMPLEX DATA EASILY

AMEDTEC ECG^{pro}

TECHNICAL DATA

Pump module	Power consumption	max. 35 W
	Dimensions	100 mm x 240 mm x 180 mm (H x L x W)
	Weight	1,1 kg
Patient module	Dimensions	30 mm x 110 mm x 135 mm (H x L x W)
	Weight	1,1 kg
Patient cable	Length	Extremity electrodes: 1,2 m (1,4 m optional available) Chest electrodes: 1,0 m (1,2 m optional available)
Electrodes	Material	The electrodes meet the biocompatibility requirements and consists of an electrode body manufactured in a latex-free thermoplastic elastomer, a superabsorbent filter and an Ag/AgCl coated electrode sensor fulfilling ANSI/AAMI recommendations
Power supply	Mains voltage / frequency	100 - 240 V / 50 or 60 Hz
	Current consumption	max. 1,1 A
	Duration	About 3,5h (corresponds approx. to 50 resting-ECG and approx. 10 stress-ECG recordings)
Battery (optional)	Charging time	Approx. 2 hours
	Length of folding arm	570 / 650 mm
Stative	Table fixation	Table width 10 - 130 mm
	Classification	Class IIa according to MDD
General	Protection class	I
	Applied part	Type CF according to IEC 60601-1, with defibrillation protection
Ambient conditions	Temperature	During operation: +10 to + 40 °C, during storage: -10 to +40 °C
	Humidity	25 - 95 %, no condensation
	Ambient pressure	700 - 1060 mbar

ACCESSORIES




QuickOne® disposable electrodes
Change electrodes after each patient
1 packing unit = 1000 pcs



Quickels® Filter electrodes
Change electrodes after 10 patients
1 packing unit = 128 pcs

AMEDTEC Medizintechnik Aue GmbH
Schneeberger Straße 5 · 08280 Aue · Germany
Phone: +49 (0)3771 59 82 70
Fax: +49 (0)3771 59 82 790
info@amedtec.de · www.amedtec.de



Manufacturer:
Cardiolex Medical AB
Sundbybergsvägen 1C
SE-171 73 Solna, Schweden

CE Certificate:
0413 ISO 13485

Errors and omissions excepted
M0685 018-0617-2_02